

## The ideal system for disinfecting **all objects and accessories** in just a few seconds!

- **Disinfect mobile phones, glasses, keys, jewellery, medical material (stethoscopes, blood pressure pads, magnifying glasses, dental equipment...)**... in just 45 seconds!
- **No chemical products**
- Tested under the **AFNOR NF T72-281** standard for airborne disinfection of surfaces
- **Eliminates at least 99.99% of micro-organisms, bacteria and viruses** (including SARS-COV-2 / Coronavirus)
- The best cost/efficiency ratio available for surface disinfection
- Adapted to all types of surfaces **without side-effects or deterioration**
- An object disinfection concept offering **confidence and reassurance to your clients and personnel**



### How UV acts

UV-C rays penetrate the heart of the DNA at 254 nanometres, the optimal wavelength for eradicating microorganisms (virus, bacteria, algae, yeast, mould...), disrupting the cell metabolism through to total destruction.

## The advantages of BIO-SCAN®

- **A certified, tried and tested process**
- Treatment time of **just 45 seconds!**
- **Practical, ergonomic and easy to use** by anyone
- **No chemical products**, no risk to health or the environment
- Economical: low purchase price and fast return on investment

### BIO-UV Group

850 avenue Louis Médard  
34400 LUNEL - FRANCE  
Tel. +33(0)4 99 133 911  
contact@bio-scan.tech  
[www.bio-uv.com](http://www.bio-uv.com)

# BIO-SCAN Cube

BIO-UV  
Group



	Number of lamps	Electric power (W)	Dimensions W x H x L (mm)	Weight (kg)	Voltage (V)	Frequency (Hz)	Power consumption (W)
BIO-SCAN® Cube	2	14	500 x 170 x 420 useful surface: 338 x 80 x 290	10	230	50/60	36

2 year  
guarantee

## Technical characteristics BIO-SCAN® Cube

- **Eliminates at least 99.99% of micro-organisms, bacteria and viruses** (including SARS-COV-2 / Coronavirus)
- Cycle time: **45 seconds!**
- Treatment of small objects (mobile phones, keys, glasses, medical devices...)
- 1.5m electric cable from the socket to the machine
- On/off switch
- Touch screen which displays the number of cycles, lamp replacement pre-alarm and alarm
- Lamp lifespan: 3000 hours (lamps specially designed for surface treatment, a BIO-UV Group patent)
- **Designed and manufactured in France**



### BIO-UV Group

850 avenue Louis Médard  
34400 LUNEL - FRANCE  
Tel. +33(0)4 99 133 911  
contact@bio-scan.tech  
www.bio-uv.com

